

ARSENIC, WHOLE BLOOD

Orderable - ASBLD

Turnaround Time: 10 days

Specimen:

6 mL K2-EDTA Royal Blue Vacutainer tube

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
0.0-50.7 nmol/L	0.0-3.8 µg/L

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here:](#)

Specimen Stability:

Ambient: 22 days
Refrigerated: 14 months
Frozen: 14 months



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis