

NIOBIMUM,
ERYTHROCYTES

Orderable - NBERC

Turnaround Time: 10 days

Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

SI Units (Reported on Patient Chart):	Conventional Units:
0.0-0.27 nmol/L	0.00-0.025 µg/L

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):

Specimen Stability:

Ambient: 22 days
Refrigerated: 14 months
Frozen: 14 months



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
ICP-MS/MS



Test Schedule:
Batched analysis