

HLA B*58:01
Allopurinol
HYPERSENSITIVITY
CASE

Orderable - HLA B*58:01 Allopurinol Hypersensitivity case

Turnaround Time: 14 days

Specimen:

2 x 4 mL Lavender (EDTA) top Vacutainer tubes

Collection Information:

Samples that are being shipped into the laboratory can be drawn Monday- Friday

Reference Ranges:

See report

Storage and Shipment:

Keep EDTA blood at room temperature and deliver within 72 hours.



Laboratory:
Transplant Lab



Requisition:
[TRANSPLANT
LABORATORY
REQUISITION](#)



Method of Analysis:
HLA Class I and Class II
 typings are performed
 using molecular
 techniques (RSSO or SSP)



Test Schedule:
As required Monday to
Friday except Stat
holidays