

ZINC
PROTOPORPHYRIN,
WHOLE BLOOD

Orderable - ZPP

Turnaround Time: 10 days

Specimen:

6 mL K2 EDTA Royal Blue Vacutainer tube

Collection Information:

[Trace Elements Collection Requirements](#)

Reference Ranges:

$\mu\text{mol/mol heme}$
0-70

Reference Ranges are based on Non-Occupationally exposed population.

Interpretive Comments:

Find Interpretive Comment and Clinical Information [here](#):



Laboratory:
Trace Elements Lab



Requisition:
TRACE ELEMENTS
REQUISITION



Method of Analysis:
Proto Fluor



Test Schedule:
Batched analysis