

MITOTANE

Orderable - MITO

Turnaround Time: 14 days

Specimen:

6 mL Red top Vacutainer tube
Avoid blood collection tubes with separator gels (e.g., Gold top)

Collection Information:

Minimum volume required is 1 mL of serum.

Reference Ranges:

7 - 90 µg/mL
Therapeutic range: 7-90 µg/mL following daily mitotane doses of 5-15 grams

Storage and Shipment:

Store and send cold or frozen.

Specimen Stability:

Ambient: 14 days
Refrigerated: 14 days
Frozen: 28 days



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Gas Chromatography (GC)



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)