

HERPES  
SIMPLEX  
VIRUS (HSV)-  
1 AND HSV-2  
IgG, SERUM

Orderable – **RESTRICTED ORDERABLES**  
(LDHSV1IGG and LDHSV2IGG)

Turnaround Time: 2 days

Alternate Name(s):

HSV-1 IgG

HSV-2 IgG

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 2 x 0.5 mL Red top Microtainers 2-10 years: 3.5 mL Gold top Vacutainer tube
Light Green (Li-Heparin) or Lavender (EDTA) top tubes are <b>NOT</b> acceptable	

Collection Information:

The minimum volume of serum required is 400 µL for adult or 300 µL for pediatric samples.

Reference Ranges:

See report

Comments:

These orderables are mainly for transplant patients. For all other patients, the HSVG orderable should be used.

Storage and Shipment:

Store and ship serum refrigerated.



**Laboratory:**  
Endocrinology Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Diasorin  
Chemiluminescence



**Test Schedule:**  
Monday - Friday  
0800-1600