

# Group and Screen (Blood Grouping, Antibody Screen and Crossmatching)

## ORDERABLE: GROUPANDSCREEN (LHSC/SJH)

**Test Details** 

Additional Information For Lab Staff Only

### TEST DETAILS

#### **ALTERNATE NAME (S):**

N/A

#### SPECIMEN:

Age	Specimen Type
Pediatric <4 months	1 mL Lavender (EDTA)
Infants >4 months and <3 years	2 mL Lavender (EDTA) top
Children >3 years and Adults	6 mL Pink (EDTA) top Vacutainer tube

#### **COLLECTION INFORMATION:**

N/A

#### **CRITICAL INFORMATION REQUIRED:**

N/A

#### TURNAROUND TIME FROM WHEN SPECIMEN IS RECEIVED IN TESTING LABORATORY:

Routine	STAT
2 hours	45 minutes

#### **REQUISITION:**

Electronic order
LABORATORY:
Blood Transfusion
TEST SCHEDULE:
As required
ADDITIONAL INFORMATION
REFERENCE RANGE:
See report
CRITICAL VALUE:
N/A
INTERPRETIVE COMMENTS:
Crossmatched blood can be available within 5 to 10 minutes upon receipt of request to convert Group and Reserve to crossmatch.
METHOD OF ANALYSIS:
Anti-Human Globulin test
FOR LAB STAFF ONLY
SPECIAL PROCESSING:
N/A
SPECIMEN STABILITY:
N/A
STORAGE AND SHIPMENT:
Transport to the Blood Transfusion Laboratory within 24 hours of collection at room temperature.
REFERRED OUT:
N/A

**COMMENTS:**