



Pathology and Laboratory Medicine

7DEHYDROCHOLESTEROL

Orderable - 7CHOL

Turn Around Time: 30 days

Specimen:

4.5 mL Light Green top Vacutainer tube

Collection Information:

Minimum volume required is 1 mL of serum.

Minimum volume required for pediatric sample is 0.2 mL.

Fasting specimen preferred. Protect specimen from light. Specimen must be labelled inside and outside light-protecting wrap.

0

Laboratory: Core Lab

Requisition: GENERAL LABORATORY REQUISITION

Ï

Method of Analysis: GC-MS

Reference Ranges:

Normal:	<5 μmol/L
Indeterminate:	5-20 μmol/L
Indicative of Smith-Lemli-	≥20 µmol/L
Opitz Syndrome:	



Test Schedule:

Referred out Monday-Thursday

Storage and Shipment:

Send sample frozen and keep protected from light at all times.



Referred Out Location: In-Common Laboratories