ALPHA2-SERUM

Orderable – A₂M

Turn Around Time: 20 days

Alternate Name(s):

A2M

Specimen:

Adult	Pediatric
5 mL Gold top	2-10 years: 3 mL Gold top
Vacutainer tube	Vacutainer tube

Collection Information:

Minimum volume of serum required is 1 mL.

Reference Ranges:

1.02 - 2.59 g/L



Test Schedule: Referred out Monday-Thursday

Storage and Shipment:

Store and ship serum frozen.

Referred Out Location: In-Common Laboratories



Pathology and Laboratory Medicine



× ||||

Laboratory: Core Lab

Requisition:

REQUISITION

Nephelometry

GENERAL LABORATORY

Method of Analysis: