

BILE ACIDS, SERUM

Orderable - BILACID

Turn Around Time: 4 days

Alternate Name(s):

Bile Acids: Total

Specimen:

5 mL Gold top Vacutainer tube

4.5 mL Green (Li-Heparin) top Vacutainer tube is also acceptable

Collection Information:

Minimum volume required for testing is 0.5 mL.

Fasting specimen preferred.

Spin, separate and store until tested.

Reference Ranges:

0.0-10.0 umol/L (fasting sample)

Storage and Shipment:

Stable 7 days at 4° C, store frozen -20° C up to 3 months



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Colourimetric Enzymatic
Assay



Test Schedule:
Batch testing - Twice per
week (Tuesday and
Thursday)