

CA125, FLUID

Orderable – CA125F

Turn Around Time: 4 hours

Specimen:

Fluid

Collection Information:

Collect a minimum of 1 mL of fluid and transfer into a Red top Vacutainer tube. Grossly hemolyzed samples are unacceptable.

May also wish to collect a Light Green top (Li-Heparin) Vacutainer tube of blood to compare fluid CA125 result to plasma CA125 result.

Reference Ranges:

No reference range available for fluid

Interpretive Comments:

It may be helpful to interpret fluid CA125 result relative to plasma CA125 result.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Roche
Electrochemiluminescence



Test Schedule:
As required