SERUM TOTAL

Orderable - CARB

Turn Around Time: 24 hours

Alternate Name(s):

Anticonvulsant Tegretol

Specimen:

Adult	Pediatric
4.5 mL Green (Lithium	0-2 years: 0.5 mL Green
Heparin) top Vacutainer	Microtainer
tube	2-10 years: 2 mL Green top tube
Serum from a 5 mL Gold top or 6 mL Red top Vacutainer	
tube is also acceptable	

Collection Information:

Collect blood aseptically in a Vacutainer tube.



Test Schedule: As required

Other times by prior arrangement.

Reference Ranges:

17 - 51 µmol/L

Critical Value: Toxic: >64 µmol/L

Comments:

CRITICAL VALUE to be phoned to Nurse or Physician immediately.



Core Lab

Requisition: GENERAL LABORATORY REQUISITION

Method of Analysis:

Kinetic Interaction of Microparticles in Solution

(KIMS)

Laboratory:



Pathology and Laboratory Medicine