

CYANIDE

## Orderable – GENX: then specify test

Turn Around Time: 25 days

---

### Specimen:

4 mL Lavender or Royal Blue top Vacutainer tube

### Collection Information:

Minimum volume required for testing is 4 mL.

Collect EDTA whole blood in contaminant-free tube.

Send in metal-free polypropylene vial.

### Reference Ranges:

< 5  $\mu\text{mol/L}$

### Special Processing:

Clearly write "CYANIDE" on the tube labels and requisition.



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Test Schedule:**  
Referred out Monday -  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)