

ESTRONE,
SERUM/
PLASMA

Orderable - ESTRO

Turnaround Time: 10 days

Alternate Name(s):

E1

Specimen:

5 mL Gold top Vacutainer tube
Red, Green (Lithium Heparin), or Lavender (EDTA) top tubes are also acceptable



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS

Collection Information:

Minimum volume of serum or plasma required is 1 mL.
Separate serum or plasma from cells within 2 hours of collection.

Reference Ranges:

Please note: New testing methodology and reference intervals as of August 30, 2022



Test Schedule:
Referred out Monday-
Thursday



Referred Out Location:
[In-Common Laboratories](#)

Male	
<= 14 days:	Not Established
Stage I (> 14 days and prepubertal):	<= 60
Stage II:	<= 82
Stage III:	37-93
Stage IV:	37-171
Stage V:	37-222
>= 19 years:	37-222

Female	
<= 14 days:	Not established
Stage I (> 14 days and prepubertal):	<= 108
Stage II:	37-122
Stage III:	55-160
Stage IV:	59-285

ESTRONE,
SERUM/
PLASMA

Female	
Stage V:	63-741
>=19 years and premenopausal:	63-741
Postmenopausal:	25-148

Special Processing:

Store and send frozen.