

ETHYL  
GLUCURONIDE

## Orderable -

Turn Around Time: 12 days

Monday/Tuesday collections received in the lab by 10:00 will have a TAT of 7 days

### Alternate Name(s):

ETG

### Specimen:

Random urine

### Collection Information:

Random urine, 5 mL minimum volume is required.

### Reference Ranges:

Ethyl Glucuronide	
Negative	<100 µg/8.8 mmol Cre
Indeterminate	100 – 200 µg/8.8 mmol Cre
Positive	>200 µg/8.8 mmol Cre
Ethanol - Urine	
≤2.4 mmol/L	

### Special Processing:

Bacterial contamination of the specimen should be prevented by collecting or transferring in sterile and air-tight containers. Store and send cold or frozen.

### Comments:

Testing is available to Liver Transplant patients and requires Biochemist Approval

### Storage and Shipment:



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION




**Method of Analysis:**  
LC-MS/MS



**Test Schedule:**  
Referred out Monday-  
Thursday



**Referred Out Location:**  
[In-Common Laboratories](#)



ETHYL  
GLUCURONIDE

Store and send frozen.