

GLIADIN
(DGP)
ANTIBODIES

Orderable - GLIAB

Turn Around Time: 5 days

Alternate Name(s):

Anti-Deamidated Gliadin Peptide

Specimen:

6 mL Red top Vacutainer tube

Collection Information:

Minimum volume of serum requested is 1 mL.

Note:

This test is not to be used for celiac disease screening unless the patient is TTGAb negative and IgA deficient. **This test is approved for Dr. Andreanne Zizzo.** All others will require biochemist approval to order this test; approval can be requested by calling the clinical biochemist at extension 77626 or the biochemist on-call.

Reference Ranges:

For Deamidated Gliadin Antibody IgA and IgG:	
<20 CU:	Negative
20-30 CU:	Weakly Positive
>30 CU:	Positive

Interpretive Comments:

The test done for Celiac Disease is Tissue Transglutamase Antibodies.

For screening for diagnosis of Celiac Disease, the patient should not be on a gluten free diet.

False negative results may occur in individuals who are IgA-deficient.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
CLIA



Test Schedule:
Referred out Monday -
Thursday



Referred Out Location:
[In-Common Laboratories](#)



GLIADIN
(DGP)
ANTIBODIES

Storage and Shipment:

Store and ship frozen.