



HIGH SENSITIVITY TROPONIN T

Orderable - TNT

Turn Around Time: 4 hours

STAT: 1 hour

Alternate Name(s):

hs-TnT
Troponin T - High sensitivity
TNT - High sensitivity

Specimen:

Adult	Pediatric
4.5 mL Green top Vacutainer tube	2-10 years: 3 mL Green top tube

Reference Ranges:

≤ 14 ng/L

Interpretive Comments:

Order hs-TnT on patients with symptoms of myocardial ischemia to diagnose AMI. Reference limit is based on 99th percentile value in healthy population. hs-TnT may be significantly elevated in non-ACS patients. Evaluate hs-TnT results in the clinical context and not in isolation.

Algorithm Details can be found

https://lhsc.omni-assistant.net/lab/Document/DocumentDownloader.aspx?Df_Guid=e595a48b-3d88-4f5d-97ac-07a64fc0fd22

Comments:

Please note that the hs-TnT reporting units are ng/L.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence



Test Schedule:
As required