



IMMUNOGLOBULIN  
D, SERUM-IGD

# Orderable - IGD

Turn Around Time: 10 days

## Alternate Name(s):

IgD

## Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Red top Microtainer 2-10 years: 2 mL Red top tube

## Collection Information:

Minimum volume of serum required is 1 mL

Lipemic or hemolysed samples are unsuitable for analysis.

Test available to LHSC/SJH physicians only.

## Reference Ranges:

7.7-132.1 mg/L

## Storage and Shipment:

Store samples at 2-8 degrees Celsius for up to 48 hours. Store specimen at -20 degrees Celsius if longer storage is required.

Intercity:

Ship frozen on ice packs



**Laboratory:**  
Core Lab



**Requisition:**  
GENERAL LABORATORY  
REQUISITION



**Method of Analysis:**  
Nephelometry



**Test Schedule:**  
Referred out Thursdays  
as required



**Referred Out Location:**  
[Laboratory Reference  
Centre](#)  
Hamilton General  
Hospital  
Core Lab - Level 1  
237 Barton Street East  
Hamilton, ON  
L8L 2X2  
Tel: 905-577-1477