

# BIOTINIDASE, PLASMA

## Orderable - BIOTP

### Specimen:

Adult	Pediatric
6 mL Green (Sodium or Lithium Heparinized) top Vacutainer tube	0-2 years: 2 x 0.5 mL Microtainers 2-10 years: 3 mL Green top tube



**Laboratory:**  
Biochemical Genetics Lab



**Requisition:**  
[BIOCHEMICAL GENETICS  
LAB REQUISITION](#)

### Collection Information:

Minimum volume of plasma required is 0.5 mL.

Store plasma at -20°C and send frozen to Biochemical Genetics.

### Reference Ranges:

6-12 umol/min/L plasma



**Test Schedule:**  
As required