

CHIKUNGUNYA VIRUS SEROLOGY

Orderable - CHIKSER

Turnaround Time: 10 days

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[Chikungunya Virus](#)

Specimen Type	Minimum Volume	Collection Kit
Serum	1mL	Gold Vacutainer tube (SST)

Critical Information Required:

Ensure the Arbovirus Information Intake form includes the following:

- Relevant travel history to endemic areas including dates of departure and return

Specify **Chikungunya** on the Arbovirus Information Intake form.

Storage and Shipment:

The specimen must be refrigerated.

The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
[Arbovirus \(Non-Zika\)
Information Intake Form](#)
This requisition **MUST**
accompany the specimen



Method of Analysis:
ELISA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Ontario
Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310