

GYNAECOLOGICAL
LIQUID
BASED PAP
FOR
CYTOLOGY

Orderable – E-order/Requisition

Turnaround Time: 1 week

Alternate Name(s):

Liquid based PAP test for Cytology
PAP test
LBP

Specimen:

Cervical, endocervical, vaginal sample in a specimen container

Collection Information:

Click here for the [External Client Cytology Collection Method](#) if you are not at LHSC or SJHC.

Reference Ranges:

See report

Interpretive Comments:

Clinical history is an important component for diagnostic interpretation.

Special Processing:

The specimen in Thinprep processed so the specimen should be collected in specific gyn. collection containers.

Critical Information Required:

Specimen container must be labelled with the patient's name.



Laboratory:
Cytopathology Lab




Requisition:
[CYTOPATHOLOGY
REQUISITION- NON-
GYNAECOLOGICAL AREA](#)



Method of Analysis:



Test Schedule:
Weekdays



GYNAECOLOGICAL
LIQUID
BASED PAP
FOR
CYTOLOGY

Unlabelled specimen containers cannot be processed.

Storage and Shipment:

Transport all specimens in a biohazard transport bag with the specimen in the zipped area, separating the requisition in the side pocket.