

HEPATITIS
DONOR
TRANSPLANT
SCREEN

Orderable - DONOR

STAT: 4 hours

Alternate Name(s):

Donor Screen

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 0.5 mL Gold top Microtainer 2-10 years: 2 mL Gold top tube

Collection Information:

DONOR careset is always a STAT request.

Reference Ranges:

See individual tests

Interpretive Comments:

Sample will be tested for Hepatitis B surface antigen, Hepatitis B core antibody (total), Hepatitis C antibody, HIV, Cytomegalovirus, RPR (screening test for syphilis), HTLV 1/2 and West Nile Virus.

Comments:

Tested only at UH Core Lab.
UH only processes TGLN/samples from deceased donors.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Electrochemiluminescence
Immunoassay



Test Schedule:
As required