

PROTEIN
ELECTROPHORESIS,
24 HOUR URINE

Orderable - PELU24

Turnaround Time: 10 days

Alternate Name(s):

PELU24
Monoclonal Protein Screen
Light Chain Screen
Bence Jones Protein Screen
Urine Protein Electrophoresis



Laboratory:
Clinical Immunology Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Sebia gel electrophoresis



Test Schedule:
Monday - Friday
0800-1600

Specimen:

24-hour urine

Collection Information:

Collection of 24-Hour Urine Sample

For 24-hour urine:
Collect urine over a 24-hour period in a sterile container with no preservatives

Please complete the laboratory label attached to the collection bottle to include patient's name, collection start time and date and collection end time and date.

Special Processing:


Minimum volume of urine required is 15 mL.

Reference Ranges:

Interpretative analysis

Interpretive Comments:

Abnormal bands are identified by immunofixation electrophoresis.



PROTEIN
ELECTROPHORESIS,
24 HOUR URINE

Comments:

Random urine specimens are no longer accepted for protein electrophoresis.

Storage and Shipment:

Store and ship urine refrigerated.