

RISPERIDONE
SERUM/
PLASMA

Orderable - RISP

Turnaround Time: 7 days

STAT:

Alternate Name(s):

Risperidol

Specimen:

6 mL Red top Vacutainer tube

or

4.5 mL Green (Lithium Heparin) top Vacutainer tube

Collection Information:

Minimum volume required is 0.5 mL.

Collect blood aseptically in a Vacutainer tube.

Reference Ranges:

Risperidone	9-Hydroxyrisperidone
19-49 nmol/L	19-211 nmol/L

Storage and Shipment:

Separate and freeze. Ship on ice packs or frozen.



Laboratory:
Toxicology/Special
Chemistry Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
LC-MS/MS



Test Schedule:
Monday--Friday
0800-1600