

TRYPsin/
TRYPsinOGEN,
SERUM

Orderable – TRYP

Turnaround Time: 2 weeks

Alternate Name(s):

Immunoreactive trypsin
Immunoreactive trypsinogen
Trypsin-like immunoreactivity

Specimen:

Adult	Pediatric
5 mL Gold top Vacutainer tube	0-2 years: 2 x 0.5 mL Red or Gold top Microtainers 2-10 years: 3 mL Red top Vacutainer tube
Light Green (Li-Heparin) or Lavender (EDTA) top tubes are NOT acceptable	

Collection Information:

Biochemist approval is required for this test. Contact the on-call biochemist prior to ordering the test. This test is not orderable in PowerChart and must be ordered with a paper requisition form.

Minimum volume of serum required is 0.3 mL.

Grossly lipemic samples will be rejected.

Reference Ranges:

180.5 – 885.3 ng/mL

Storage and Shipment:

Store and ship serum frozen.



Laboratory:
Core Lab



Requisition:
GENERAL LABORATORY
REQUISITION



Method of Analysis:
Radioimmunoassay



Test Schedule:
Referred out Monday -
Friday



Referred Out Location:
[LifeLabs](#)