

LYMPHOCYTE
PROLIFERATION

Orderable - LYMPRO

Alternate Name(s):

Lymphocyte Mitogen

Specimen:

2 X 6 mL Dark Green top Vacutainer tubes (Sodium Heparin)

Collection Information:

A control specimen must also be sent. Control specimen should be labeled "Control Sample for Patient _____".

Patient and control samples must be sent together.

Pricing: \$490 per test plus \$38 for shipping.

This test is only available at the request of Dr. Ashley Geerlinks, Dr. Mark Kuprowski, Dr. Harold Kim, Dr. Samira Jeimy and Dr. Rongbo Zhu. If any other physicians request this test, please contact the Biochemist on call for authorization. Any other physicians ordering this test need a full immunology consult to determine if it is necessary.

Reference Ranges:

See report

Special Processing:

Sample must arrive at the testing lab by 10:00 am on the day of the test (Thursday), no greater than 24 hours after the initial blood draw, accompanied by a control blood sample and a completed requisition. Please notify the testing lab before sending samples.

Specimen should be drawn first thing in the morning before 10:00 and shipped the same day. Sample must be less than 24 hours old to be tested. A control specimen must accompany the patient specimen.



Laboratory:
Core Lab



Requisition:
[SICK KIDS LYMPHOCYTE
PROLIFERATION - PHA
REQUISITION](#)



Test Schedule:
Sample must arrive at the testing lab on Thursday by 10:00 with a completed requisition.



Referred Out Location:
Hospital for Sick Children
Division of Immunology
and Allergy - Service Lab

Attention: Dr. Harjit Dadi
Peter Gilgan Centre for
Research and Learning
686 Bay Street
Room: 219430
Toronto ON M5G 0A4
Tel: 416-813-6254



LYMPHOCYTE PROLIFERATION

Be sure that shipping is aware that the package must be delivered by 0900 the following morning. This requires the package to be shipped via Purolator and there is an additional \$38.00 cost; but the sample must be shipped for 0900 delivery. Please ensure you get a tracking number.

Critical Information Required:

Special requisition from Hospital for Sick Children must be completed in full and accompany specimen.

Storage and Shipment:

Ship whole blood at room temperature.