

SYPHILIS SEROLOGY

Orderable - SYPHILIS

Turnaround Time: Up to 6 days

Specimen Type and Collection Information:

Please review the Public Health Ontario Laboratory Test Information Index:

[Syphilis - Serology](#)

If ordering as part of a Prenatal Serology work-up, see Prenatal Serology.

Specimen Type	Minimum Volume	Collection Kit
Serum	2.0 mL	Blood, clotted – Gold Vacutainer tube (SST)

For initial syphilis testing on newborns, it is recommended that maternal serum also be submitted with a separate requisition.

Storage and Shipment:

The specimen must be refrigerated.
The specimen must be received within 24 hours of collection.



Laboratory:
Microbiology



Requisition:
See electronic order



Method of Analysis:
CMIA
RPR
TPPA



Test Schedule:
Referred out weekdays



Referred to:
Public Health Laboratory
1200 Commissioners Rd E
London, Ontario
519-455-9310